

TM

ergonomic productsTM
Performance Without Compromise[®]



SlimLine Universal HeadwallTM

The following pages excerpted from Ergonomic Products catalog: 950-075revM

[Click here to download full catalog](#)

www.ergonomic-products.com

SlimLine Universal Headwall™ (Item C3A1)



Monitor and Slat Wall accessories not included.

Our new array of **SlimLine** products gives you **all the EP productivity** you've come to expect – in a **streamlined, modern styling**.

Additionally, the shallower Corian Work Surface on the SlimLine Universal Headwall makes it ideal for operatories with restricted room-length, while still providing ample counter space.

Standard Features

Upper Storage Cabinet

- 3.7 cubic feet of storage
- Full-swing door on upper Computer Cabinet, opens to Adjustable Shelves.
- Cabinet Shelf (below Computer Cabinet) includes double Slide-Doors for gloves and other bulk storage (see below for optional configurations)
- Monitor Mounting Bracket (monitor not included)
- Built in LED under-cabinet Lights

Integrated Slat Wall

- Nearly 4 sq ft of Slat Wall to display and store accessory materials (see option below)

Work Surface and Utilities

- 2.3 sq ft of solid Corian Work Surface
- Open Shelf suitable for Zirc Tub placement
- Quick Connect Umbilical port (standard)
- Utility Access Panel

CONFIGURATION and UPGRADE CHOICES

Cabinet Color/Finish Selection:

- ☐ Cabinet available in standard EP color choices.
- ☐ - or - Upgrade to woodgrain finish. (Standard EP selections)

Corian Work Surface:

- ☐ Available in standard and premium EP Corian color choices.

Utility Configuration:

- ☐ Standard Umbilical Quick Connect
- ☐ - or - Utility Plate and Clamp-mount for 3rd party utilities

Cabinet Shelf Configuration Options (shown right):

- ☐ Open Cabinet Shelf, removes Slide-Doors for easy access
- ☐ - or - No Cabinet Shelf, increases available Slat Wall space by nearly 2 additional sq ft



Open Cabinet Shelf




No Cabinet Shelf, for Expanded Slat Wall





Dimensions & Utilities



⚡ Electrical Specifications

120v AC, 60HZ, 8A
Cat 5 (or better) Data Jack

🌀 Air / Vacuum Specifications

Air: 1/4" Dental Style Female QD, 80 psi
Vacuum: 1/2" Copper or PVC Pipe

💧 Plumbing Specifications

Dental Line feed:
- 1/2" Cold or Dedicated Dental supply line with straight 1/4 turn-stop and 3/8" compression

Ergonomic Products offers a range of **Slat Wall** items including:

- ☐ Med-waste Bracket
- ☐ Remote Control Holder
- ☐ Sanitation Wipes Holder
- ☐ Pad & Pen Holder
- ☐ "Goodie Bag" Brackets
- ☐ Printer Shelf
- ☐ Monitor Bracket
- ☐ Glove Box Holders

See page 116 for details.

Ideal Room Location

Similar in function to our standard headwall, the SlimLine Universal Headwall handles storage and material presentation with great efficiency and is most effective when used in a dual / head-entry format.

However, the streamlined Corian Work Surface allows the unit to be located closer to the patient chair - allowing for a shorter overall room length.

