

TM

# ergonomic products<sup>TM</sup>

Performance Without Compromise<sup>®</sup>



## RapidCarts<sup>TM</sup>

The following pages excerpted from Ergonomic Products catalog: 950-075revM

[Click here to download full catalog](#)

[www.ergonomic-products.com](http://www.ergonomic-products.com)


# RapidCart™ Systems

Don't duplicate expensive equipment – **mobilize it!**

RapidCarts allow instantaneous deployment and setup of less frequently used treatment instruments such as those in endodontics, surgery and advanced diagnosis. Distributing technology from room to room – rather than targeting treatment to individual rooms – creates a more flexible and responsive clinical environments and is one of Ergonomic Products' core efficiency principles:

## Any Procedure, Any Room, Any Time.

Ergonomic Products' line of RapidCarts has been designed specifically for dental practices with a narrow profile for ease of movement in operatories, and with large casters for added stability. The standard product Work Surface is durable Corian, to stand up to the rigors of a medical environment.



Mobile Cart Bay in a 7 Chair Practice

## COLOR OPTIONS - ALL CARTS

### Work Surface:

- ☐ EP Standard Color Corian
- ☐ - or - Premium Corian upgrade
- ☐ - or - Powder-coated top (reduced price)

### Cabinet Color Selections:


- ☐ Exterior finish available in standard EP colors
- ☐ - or - upgrade to woodgrain finish  
(Standard EP woodgrain finishes)

## Mobile Storage RapidCart™ (item R1C1)



### Standard Features

- Anodized aluminum Corner-framing
- 4 milled Handpiece Holders
- High-impact plastic Casters
- Precision-bearing Drawer Slides
- Soft-touch Drawer Close for secure transport



Height, width and depth dimensions above are standard for the three carts shown here.

See next page for shelving and utility configurations.





## Doctor Specialty RapidCart™ (item R1A1)



### Standard Features

- Wiring Apertures on Work Surface, with stainless steel Cover Plates
- 5 Outlet Power Strip mounted on rear
- Velcro cord management Loop
- 2 Side Aperture Wiring Openings
- Adjustable Shelf
- Milled Small Tool Recesses on Work Surface.
- Two deep Drawers
- Soft-touch Drawer Close for secure transport

### Suggested Accessories (see pages 114 & 117)

- 1 Zirc Tub for Shelf (item M1E1 or M1E2)
- 2 Lewis Bins, in Drawers (item M1D1)
- Foot Pedal Holder (2 sizes)

## Rapid Deployment RapidCart™ (Item R1B1)



### Standard Features

- Fixed Slide Brackets to receive Standard Zirc or other supply tubs
- Two deep Drawers
- Soft-touch Drawer Close for secure transport

### Suggested Accessories (see page 114)

- Standard Zirc Tub (item M1E2)
- 2 Lewis Bins, in drawers (item M1D1)