



TM

ergonomic productsTM
Performance Without Compromise[®]



Patient & Clinician Seating

The following pages excerpted from Ergonomic Products catalog: 950-075revM

[Click here to download full catalog](#)

www.ergonomic-products.com

3-Series™ Patient Chair

Heat, Massage & Traverse for the best possible patient and clinician experience.

When your patients' comfort is job one, why would you want anything less than the **most comfortable, feature-rich** chair on the market?



CONFIGURATION

Function:

- ☐ EP 390 - Heat, Massage & TraverseItem S1A1
- ☐ EP 370 - Traverse OnlyItem S1A2
- ☐ EP 330 - Heat & Massage OnlyItem S1A3
- ☐ EP 310 - No Heat, Massage or TraverseItem S1A4

OPTIONS

- ☐ Self-tensioning Arm Slings
- ☐ 2nd Chair Control Touchpad
- ☐ Foot Control Pedal
- ☐ 2" Accessory Mounting Post
- ☐ Junction Box

COLOR CHOICES & STITCHING

Upholstery:

- ☐ Seamless, for a sleek, easy-to-clean appearance well suited to the dental environment (as shown above)
- ☐ - or - Upgrade to Decorative Contour Stitching as shown on following page (Tri-tone color not currently available)

Upholstery is available in over 60 shades of high-performance Ultraleather. These are our some of our most popular choices, but we have dozens more to suit any pallet!

Swatch-cards available.



ergonomic products
Performance Without Compromise™




www.ergonomic-products.com • 866-374-6487

Excerpted from EP Catalog:

950-075revM



Dimensions, Features & Specifications



Absolutely the most feature-rich chair on the market, the EP 390 was designed to consolidate the best-in-class features appearing over the past 3 decades.

With focused attention to the details of patient comfort and practitioner ease-of-use, the EP 390 will redefine your perspective on seating.


Standard Features

- Heated base
- Lower lumbar massage
- 10-1/4" Traverse length for perfect positioning
- Memory Foam support
- Rich Ultraleather® upholstery
- Hydraulic lift
- Audio inputs

Electrical Specifications

120v AC, 6 amp fuse

Audio line (if desired)



Vision Patient Chair (Item S1B2)

Every practice deserves a chair that puts **patient comfort first**, without pushing past the limits of your **budget**.

That's why we created the **Vision Patient Chair**.

The Vision offers all the essential functions of any mainline chair—while being more comfortable, requiring less maintenance, and coming in at a very attractive price.

Built to Ergonomic Products exacting production standards, the Vision Patient Chair is a durable, high-quality unit that will provide years of exceptional service in your practice.




OPTIONS


☐ Narrow Back version (see opposite page)
with included Self-tensioning Arm Slings

☐ 2nd Chair Control Touchpad
☐ Wireless Foot Control Pedal


COLOR CHOICES




Jet



Gray



Parchment



Naval



Tannery

Other colors may be available.
Call for details.



ergonomic products
Performance Without Compromise™

www.ergonomic-products.com • 866-374-6487

Excerpted from EP Catalog:

950-075revM



Dimensions, Features & Specifications



Narrow Back
option with included
Self-tensioning
Arm Slings.

The Vision Patient Chair is a remarkably versatile entry-level seat featuring the same build quality of our premium chair.

With the added comfort of our patented Self-tensioning Arm Slings (available on the Narrow Back configuration), this chair provides a patient experience that far exceeds its modest price.

Standard Features


- Powered Backrest and elevation
- Multi-point articulating Headrest
- Programmable Control Touchpad at arm
- 2 presets plus exit position
- Easy locking swing-up Arms
- Durable 2.5" thick, extra firm foam

Electrical Specifications

120v AC, 3 amp fuse



Shown with optional
second Chair Control
Touchpad



SpriteTx™ Pediatric Lift Bench (Item S5A1)



*Shown with optional
A/W Syringe & Staging Tray*

Versatile, comfortable
delivery that will **elevate**
your patients...and your
productivity!

Our signature Ergonomic Products **compact delivery channel** and an **adjustable elevation range of 16"** allows much closer positioning to the patient than other benches.

The swiveling channel and assistant arm are equally adept at both four- and two-handed delivery – *and* allow for instant **left/right convertibility**.

An optional **staging tray** provides additional work surface for more extensive treatment, but folds down and away when not in use.

And the rich, durable **Ultraleather®** upholstery is available in a range of bright and **appealing colors** to wow your patients!

Standard Features

- 1 High-speed Air
- 1 Slow-speed Air
- 1 HVE
- 1 Saliva Ejector
- 1 Air/Water Syringe
- 2 Liter Self-contained Water Supply
- Foot Pedal control
- Removable Head & Neck Support
- Rich Ultraleather upholstery
- Dual-density Foam for better and more secure patient support
- Durable Mechanical Lift with igus® Carrier System


COLOR CHOICES

Upholstery:


- ☐ Available in fun, kid-friendly Ultraleather colors
- ☐ - or - Call for additional options. Swatch cards available.



291-2527 Lagoon



291-4460 Parrot



291-5250 Sassy



291-2685 Bayou

UPGRADES

Utilities:

- ☐ 2nd HVE
- ☐ 2nd Air/Water Syringe on Doctor's Channel
- or - 3rd Handpiece (High or Slow-speed)
- ☐ Staging Tray on Swivel-Arm

Fiber Optics:

- ☐ Fiber Optics Light Source & Power Pack


Configuration:

- ☐ Removable magnetic Staging Tray on swivel arm





Dimensions, Features & Specifications



Utility Connections



Integrated floor box for power, air and vacuum connections

Electrical Specifications

120v AC, 3 amp fuse

Air / Vacuum Specifications

Air: 1/4" Compression, Female, 80 psi
Vacuum: 1/2" Copper or PVC stubbed

Plumbing Specifications

N/A - Integral Bottle only



Clinician Seating - RGP


Operator seating is often the most neglected piece of dental technology – yet in many ways it's the **most important!**

Today, with the wide range of services offered, many practitioners find themselves in a fixed seated position for extended time periods — a situation that is extremely deleterious to the spine and to physical health in general.



Finding the right seat will contribute immeasurably to your long-term health.

Watch Dr. Ahearn's [short video](#) on posture, seating options, what makes a great choice in clinical seating.



RGP 400-D LINE

The patented design of the **RGP 400** line helps protect your most important asset:
your health.

Standard Features

- Patented split Backrest for lower lumbar support
- Contour Seat Cushion
- Built with standard cylinder height range (18.5"– 24.5" seat level)
- "Active Seating" responds to weight and position changes
- Rich Ultraleather® upholstery
- Single lever, multi-function Control
- Assistant Stool includes adjustable Foot-ring for added comfort and utility



Dr./Hyg. Stool (item S2A1)



Asst. Stool (item S2A2)

Unique to the 400 series is its hydraulic mechanism which allows the stool's back and seat to 'float', following one's movement. This "Active Seating" is further enhanced by the adjustable heights of both the backrest and the stool itself.

Early in our research of the most ergonomically efficient dental methods, we realized that belly bars are not typically necessary. In fact, they often get in the way. With proper positioning of materials and supplies, assistants are able to maintain a sound ergonomic posture that provides better comfort and improved productivity. However, belly bars are available.

COLOR CHOICES & CONFIGURATION

RGP stools are available in any of our EP standard color selections to match or compliment the 3-Series Patient Chair. Swatch samples and custom colors (upgrade) are available. Please call for details.

The Doctor/Hygienist stool is available with optional armrests (shown above).



ergonomic products
Performance Without Compromise™

www.ergonomic-products.com • 866-374-6487

Excerpted from EP Catalog:

950-075revM

Clinician Seating - Crown

CROWN SEATING

For nearly 20 years **Crown Seating** has been an industry leader in the design and manufacturing of ergonomically correct Dental and Medical stools. Utilizing the **highest quality materials and craftsmanship**, Crown Seating manufactures a wide variety of styles, at **many different price levels**, for professionals across multiple industries.

Aurora



Dr./Hyg. Stool (item S2B1)



Asst. Stool (item S2B2)

Standard Features

- Triple lever, independent adjustment
- Advanced molded Contour with high-resistance molded foam cushioning
- Stable 5-leg chrome Base
- Vinyl upholstery
- Multi-angle pelvic positioning to reduce lower back strain
- Assistant Stool includes adjustable Foot-ring for added comfort and utility

Denver



Dr./Hyg. Stool (item S2B3)



Asst. Stool (item S2B4)

Standard Features

- Dual lever adjustment
- Ergo Angle style Saddle
- Stable 5-leg powder-coated Base
- Vinyl upholstery
- Multi-angle pelvic positioning to reduce lower back strain
- Seat Groove reduces tailbone pressure
- Assistant Stool includes adjustable Foot-ring for added comfort and utility

Sterling



Dr./Hyg. Stool (item S2B5)

Shown
with op
Belly B:



Asst. Stool (item S2B6)

Standard Features

- Triple lever, independent adjustment
- Deeply contoured Seat with significant "waterfall" front
- Stable 5-leg powder-coated Base
- Vinyl upholstery
- Multi-angle pelvic positioning to reduce lower back strain
- Assistant Stool includes adjustable Foot-ring for added comfort and utility

