

SpriteTx™ Pediatric Lift Bench (Item S5A1)



*Shown with optional
A/W Syringe & Staging Tray*

Versatile, comfortable
delivery that will **elevate**
your patients...and your
productivity!

Our signature Ergonomic Products **compact delivery channel** and an **adjustable elevation range of 16"** allows much closer positioning to the patient than other benches.

The swiveling channel and assistant arm are equally adept at both four- and two-handed delivery – *and* allow for instant **left/right convertibility**.

An optional **staging tray** provides additional work surface for more extensive treatment, but folds down and away when not in use.

And the rich, durable **Ultraleather®** upholstery is available in a range of bright and **appealing colors** to wow your patients!

Standard Features

- 1 High-speed Air
- 1 Slow-speed Air
- 1 HVE
- 1 Saliva Ejector
- 1 Air/Water Syringe
- 2 Liter Self-contained Water Supply
- Foot Pedal control
- Removable Head & Neck Support
- Rich Ultraleather upholstery
- Dual-density Foam for better and more secure patient support
- Durable Mechanical Lift with igus® Carrier System


COLOR CHOICES

Upholstery:

- ☐ Available in fun, kid-friendly Ultraleather colors
- ☐ - or - Call for additional options. Swatch cards available.



291-2527 Lagoon



291-4460 Parrot



291-5250 Sassy



291-2685 Bayou

UPGRADES

Utilities:

- ☐ 2nd HVE
- ☐ 2nd Air/Water Syringe on Doctor's Channel
- or - 3rd Handpiece (High or Slow-speed)
- ☐ Staging Tray on Swivel-Arm

Fiber Optics:

- ☐ Fiber Optics Light Source & Power Pack


Configuration:

- ☐ Removable magnetic Staging Tray on swivel arm





Dimensions, Features & Specifications



Utility Connections



Integrated floor box for power, air and vacuum connections

Electrical Specifications

120v AC, 3 amp fuse

Air / Vacuum Specifications

Air: 1/4" Compression, Female, 80 psi
Vacuum: 1/2" Copper or PVC stubbed

Plumbing Specifications

N/A - Integral Bottle only

