




ergonomic products™  
Performance Without Compromise®



# HYGIENE ROOM PACKAGE






## The New Standard in Hygiene Production



The compact InWall all-in-one delivery platform at the core of our Hygiene Package provides the same enhanced production found in all our Workstations.

With our EP 330 Patient Chair (featuring heat and massage), your Hygiene rooms will be attractive, comfortable spaces that will guarantee your patients return time and time again.

- High-performance, articulating delivery platform
- Space Saving InWall construction
- Material consolidation for improved ergonomics



## Performance without Compromise.®

**We all know Hygiene drives a practice;** it's the portal through which most new patients enter, and the vehicle that keeps them long-term. **To sustain growth, a practice CANNOT limit capacity in this crucial area.**

Is your hygiene set-up equal to the demands placed on it? **Ours is – and then some!**

Our 32" Hygiene InWall delivery platform features an articulating, 3-point swivel for left/right convertibility and reduced range-of-motion which is the hallmark of the EP Workstations. The concealed hydraulic pull-downs in our Resupply allow for exceptional storage capacity with minimal visual impact. Space-saving InWall Storage and Comfort accommodate both the patient experience and the demands of the busy clinician. The EP 330 Patient Chair brings the package together and is so comfortable your patients will ask where they can get one of their own!

Hygiene has a huge impact on the health of your practice; but it doesn't have to have an equally large impact on your Floorplan. **With the EP Hygiene Room, you don't need to sacrifice Productivity when space is at a premium.**

23\_01\_06





## 32" Hygiene InWall™ Workstation (items D2A1 & D2A4)

### Standard Features Include:

- 3-point Swivel for 2-hand/4-hand convertibility
- 1 High Speed Air
- 1 Low Speed Air
- 1 A/W Syringe
- 1 HVE
- 1 Saliva Ejector
- Parkell Manual Scaler - or - Parkell Turbo Sensor Scaler
- 1 EP Consumable Bin
- 1 EP Med-waste Bin
- 8 Milled Hand-piece Holders
- 2-Liter Self-contained water supply
- Auxiliary Air/Water Ports
- Bib/Glove/Product storage (D2A4)

*Open storage is recommended in Central Resupply, but is less desired in patient spaces. The Hygiene Resupply eliminates the under-storage of our standard unit and features doors for a patient-friendly appearance.*



## Hygiene Resupply Pulldown (item C1A2)

### Standard Features Include:

- Hydraulic pull-down shelving provides easy access to materials.



## 16" InWall Comfort Station (item N1A1)

### Standard Features Include:

- Large, hinged mirror covers InWall storage
- Impact resistant, seamless face-frame
- Chrome faucet
- Brushed stainless under-mount sink
- Beveled Corian counter
- Corian backsplash
- Utility access panel
- Recessed paper towel and cup dispenser behind mirror



## 16" InWall Storage (item N1B1)

### Standard Features Include:

- Impact resistant, seamless face-frame
- Bulk storage specifically designed to hold standard glove boxes

### Optional Features Include:

- Flip-up writing surface, or
- Med-waste flip-out
- Chart holder, or
- Paper-towel holder



## FastTrack LS Track Light (item L2A1)

### Standard Features Include:

- All LED lighting
- 7000-32000 lux intra-oral light
- Built-in touchless On/Off control
- Tune-able wave-length allows intra-oral use with a range of curable composites
- 360° light-arm rotation
- Single-piece aluminum extrusion body
- Milled billet aluminum trolley
- Igus® Energy Chain for reliable cable management



## EP 330 Patient Chair (item S1A3) - Standard Features Include:

- Heated base
- Lower Lumbar Massage
- Memory Foam support
- Rich Ultraleather® upholstery
- Hydraulic lift
- Audio inputs

Product shown with optional self-tensioning arm slings, and decorative contour stitching. (Tri-tone upholstery scheme not currently available.)



## RPG 400-D (item S2A1) Doctor/Hygiene Stool

- Patented split backrest for lower lumbar support
- "Active Seating" responds to weight and position changes
- Rich Ultraleather® upholstery
- Single lever, multi-function control

## Hygiene Room Package - Standard Pricing



Item #	Description	Price
<b>Delivery</b>		
D2A1	32" InWall Workstation (Lower Unit, Std Config) .....	\$ 6,499.00
D2A4	32" InWall Workstation (Upper Storage Unit) .....	799.00
<b>Storage</b>		
C1A2	Hygiene Resupply PullDown .....	1,849.00
N1B1	16" InWall Storage.....	739.00
N1A1	16" InWall Comfort Station.....	1,859.00
<b>Lighting</b>		
L2A1	FastTrack Track Light .....	3,899.00
<b>Patient &amp; Practitioner Seating</b>		
S1A3	EP 330 Patient Chair.....	9,009.00
S2A1	RGP 400-D Doctor/Hygiene Stool.....	1,129.00
<b>TOTAL</b>		<b>\$ 25,782.00</b>

Pricing does not include X-rays, CPU, or other supplies shown. All pricing based on standard configurations, with cabinetry & Corian in white. Prices subject to change. Call for details.

23\_01\_06